

X17B2 SAFETY CHECKLIST			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
Rev: A	Effective: 6-9-2003	Page 1 of 1	Number: LS-SCL-0038
Reviewed by: K. Warburton		Reviewed by: N/A	Approved by: A. Ackerman

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

Completed by: _____

Date: _____

OPCO Approval: _____

Date: _____

- ☐ 1. X17B1 Safety Checklist complete and in place.
- ☐ 2. X17B2 hutch – threshold door piece (painted yellow) bolted in place under each of 2 doors going into hutch.
- ☐ 3. X17B2 hutch – lead shielding mat against threshold door piece in each of 2 doors going into hutch.
- ☐ 4. X17B2 hutch interior – water cooled backstop/ Bremsstrahlung shield #4 downstream end of hutch (see photo)
- ☐ banded and in the beam path (**for X17B2 operation**)
- OR**
- ☐ banded and removed from beam path (**for X17B3 operation**)